

ECLIPSE

PRP[®] HC
HIGH CONCENTRATION



The New Eclipse PRP[®] 44mL Kit

The Eclipse PRP 44mL kit offers practitioners the opportunity to obtain more PRP from a single blood draw while continuing to use the industry-leading and familiar Eclipse PRP single-spin system. With minimal changes to the current smaller volume draw procedures, both staff and patients will not encounter any disruption in their normal experience.

LIST OF COMPONENTS

- (2) 22mL PRP tubes
- (1) Blood Collection Device
- (1) Latex Free Tourniquet
- (4) 10mL Syringes
- (2) Hypodermic Needle
- (2) Venting Needle
- (2) Aspiration Needle
- (2) Female to Female Luer Adapter

PRODUCT SPECS

- Eclipse PRP[®] HC 44mL Kit**
- PLT Concentration* ~4.5x
- Total Platelets ~10 x billion
- RBC <1%
- WBC <5%
- Mononuclear cells 86.2%
- PLT Yield >80% (STD Dev. +/- 9%)
- pH 7.2-7.7
- PRP Volume ~8mL

ECLIPSE[®]